

# Protocol: Quantifying species abundance using 16S Sanger Sequencing and CASEU

Jack Reddan, Kapil Amarnath

May 25, 2022

## Contents

<b>1</b>	<b>Culture Collection</b>	<b>3</b>
1.1	Overview . . . . .	3
1.2	Materials . . . . .	3
1.3	Method . . . . .	3
<b>2</b>	<b>gDNA Isolation</b>	<b>3</b>
2.1	Overview . . . . .	3
2.2	Materials . . . . .	3
2.3	Method . . . . .	3
<b>3</b>	<b>16S rDNA Amplification</b>	<b>4</b>
3.1	Overview . . . . .	4
3.2	Materials . . . . .	4
3.3	Method . . . . .	4
<b>4</b>	<b>16S rDNA Purification</b>	<b>5</b>
4.1	Overview . . . . .	5
4.2	Materials . . . . .	5
4.3	Method . . . . .	5
<b>5</b>	<b>Sanger Sequencing</b>	<b>6</b>
5.1	Overview . . . . .	6
5.2	Materials . . . . .	6
5.3	Method . . . . .	6
<b>6</b>	<b>Species Abundances with CASEU</b>	<b>6</b>
6.1	Overview . . . . .	6
6.2	Materials . . . . .	6
6.3	Method . . . . .	7

<b>A</b>	<b>PCR</b>	<b>7</b>
A.1	Master Mix Recipe . . . . .	7
A.2	Thermocycler Program . . . . .	7
<b>B</b>	<b>Genewiz</b>	<b>7</b>
B.1	Sample Mix Recipe . . . . .	7

# 1 Culture Collection

## 1.1 Overview

This section describes collecting cultures for which you wish to quantify species abundances. In addition to any experimental mixtures, pure cultures of community members are needed for the analysis.

## 1.2 Materials

- Pure and mixed liquid cultures
- 2-mL microcentrifuge tubes

## 1.3 Method

1. In a 2-mL microcentrifuge tube, Collect 1 ODmL of liquid culture for all mixtures and pure cultures. Spin-down and remove supernatant.
2. Either store pellets at  $-20^{\circ}\text{C}$ , or continue to [gDNA Isolation](#).

# 2 gDNA Isolation

## 2.1 Overview

This section describes how to isolate gDNA to be used for 16S rDNA amplification. These steps are for the QIAGEN DNeasy Blood & Tissue Kit. If you happen use a different kit, please follow its instructions, and then proceed to [16S rDNA Amplification](#).

## 2.2 Materials

- QIAGEN DNeasy Blood & Tissue kit
- 100% ethanol
- Cell pellets from cultures collected in [Culture Collection](#)

## 2.3 Method

Follow the subsequent instruction in parallel for all samples to be processed.

1. Resuspend pellet in 180  $\mu\text{L}$  of buffer ATL, and add 20  $\mu\text{L}$  of Proteinase K. Mix thoroughly by vortexing. Incubate at  $56^{\circ}\text{C}$  for 1 hour, and vortex every 15 minutes.
2. Vortex for 15s, then add 200  $\mu\text{L}$  of buffer AL and mix thoroughly by vortexing. Immediately add 200  $\mu\text{L}$  of 100% ethanol and, again, mix thoroughly by vortexing.

3. Pipet the mixture from step 2 into a DNeasy Mini spin column placed in a 2-mL collection tube. Centrifuge at  $\geq 6,000 \times g$  for 1 min. Afterwards, discard flow-through and collection tube.
4. Place DNeasy Mini spin column into a new 2-mL collection tube and add 500  $\mu\text{L}$  of buffer AW1. Centrifuge at  $\geq 6,000 \times g$  for 1 min. Afterwards, discard flow-through and collection tube.
5. Place DNeasy Mini spin column into a new 2-mL collection tube and add 500  $\mu\text{L}$  of buffer AW2. Centrifuge at  $20,000 \times g$  for 3 min. Afterwards, discard flow-through and collection tube.
6. Place DNeasy Mini spin column into a clean 1.5-mL microcentrifuge tube. Add 200  $\mu\text{L}$  of buffer AE and incubate at room temperature for 1 min. Centrifuge at  $\geq 6,000 \times g$  for 1 min to elute.
7. Check gDNA concentration on the nano spectrophotometer. Then store at  $-20^\circ\text{C}$  or continue to **16S rDNA Amplification**.

## 3 16S rDNA Amplification

### 3.1 Overview

This section describes how to amplify the 16S rDNA using Phusion polymerase for Sanger sequencing. If you are using a different polymerase, please follow its specific instructions and proceed to **Sanger Sequencing**.

### 3.2 Materials

- Phusion polymerase
- 5X High Fidelity (HF) reaction buffer
- dNTPs
- 27F forward and 1492R reverse primers
- gDNA template from **gDNA Isolation**

### 3.3 Method

Follow the subsequent instruction in parallel for all gDNA samples collected in **gDNA Isolation**.

1. Mix the ddH<sub>2</sub>O, dNTPs, forward/reverse primers, 5X HF buffer, and polymerase according to **Appendix A** for  $N + 1$  reactions, where  $N$  is the number of gDNA samples to be processed.

2. Pipet 38  $\mu\text{L}$  of the PCR master mix into 200- $\mu\text{L}$  PCR tubes. Then add 2  $\mu\text{L}$  of template gDNA, at a concentration of 100 ng/ $\mu\text{L}$ , to the mix for each sample.
3. Place reactions into the thermocycler and run the PCR program specified in [Appendix A](#).
4. Once the run is complete, check for synthesis of the desired product using gel electrophoresis. If the proper product was amplified, either store the reaction at  $-20^{\circ}\text{C}$  or proceed to [16S rDNA Purification](#)

## 4 16S rDNA Purification

### 4.1 Overview

This section describes how to purify the 16S rDNA amplicons from [16S rDNA Amplification](#) using the QIAGEN QIAquick PCR Purification Kit. If you are using a different purification kit, please follow its instructions and continue to [Sanger Sequencing](#).

### 4.2 Materials

- 16S rDNA PCR amplicon samples
- QIAquick PCR Purification Kit

### 4.3 Method

Follow the subsequent instructions in parallel for all PCR amplicons synthesized in [16S rDNA Amplification](#).

1. Add 5 volumes of buffer PB to 1 volume of PCR sample and mix by pipet.
2. Place a QIAquick spin column in a provided 2-mL collection tube. Then apply the sample to the column and centrifuge at  $\geq 17,900 \times g$  for 60 s. Discard the flow-through.
3. Add 750  $\mu\text{L}$  of buffer PE to the column and centrifuge at  $\geq 17,900 \times g$  for 60 s. Discard the flow-through.
4. Centrifuge the column at  $\geq 17,900 \times g$  for 60 s, then transfer the column to a clean 1.5-mL microcentrifuge tube.
5. To elute, add 30  $\mu\text{L}$  of buffer EB to the center of the column. Let stand for 1 min 60 s, and then centrifuge at  $\geq 17,900 \times g$  for 60 s.
6. Check the concentration of the purified product on the nano spectrophotometer. Then either store at  $-20^{\circ}\text{C}$  or proceed to [Sanger Sequencing](#).

## 5 Sanger Sequencing

### 5.1 Overview

This section describes how to prepare the purified 16S rDNA samples obtained from [16S rDNA Purification](#) for Sanger sequencing by [Genewiz](#). If you are using a different sequencing platform or provider, please follow the relevant instructions, as they most certainly will differ from what is listed here, and continue to [Species Abundances with CASEU](#)

### 5.2 Materials

- Purified PCR 16S rDNA amplicons
- 27F forward primer
- 8-tube PCR strip

### 5.3 Method

Follow the subsequent instructions in parallel for all PCR amplicons purified in [16S rDNA Purification](#).

1. Dilute an aliquot of the 16S rDNA amplicon to a concentration of 3 ng/ $\mu$ L.
2. Add the amplicon dilution, 27F forward primer, and ddH<sub>2</sub>O to a PCR tube following the recipe in [Appendix B](#).
3. Label tubes using the tube number with the prefix *HL* (e.g. HL1, HL2, ...HLN), noting what each tube contains.
4. Submit a Genewiz Sanger sequencing order online and print the order form. Place both the PCR tubes and order form in a clear, sealable plastic bag and deposit in the Genewiz drop-box.

## 6 Species Abundances with CASEU

### 6.1 Overview

This section describes how to quantify species abundances with the R package **CASEU** (Community Analysis via Sanger Electropherogram Unmixing) [1] using the sequencing results obtained from [Sanger Sequencing](#). If you choose to use different software to quantify abundances, I do not know why you are reading this protocol.

### 6.2 Materials

- .ab1 sequencing files for pure and mixed cultures
- Computer with an installation of [R](#)

## 6.3 Method

1. Refer to the vignette found on the [CASEU Repository](#) for instructions on installation and a full introduction to CASEU. Also, you can see [Jack's CASEU Repository](#) for an example script.

## A PCR

### A.1 Master Mix Recipe

### A.2 Thermocycler Program

## B Genewiz

### B.1 Sample Mix Recipe

## References

- [1] N. Cermak, M. S. Datta, and A. Conwill, “Rapid, inexpensive measurement of synthetic bacterial community composition by sanger sequencing of amplicon mixtures,” *iScience*, vol. 23, p. 100915, mar 2020.